




LBY-XC40B/ LBY-XC20B

Auto Dynamic Erythrocyte Sedimentation Rate (ESR) Monitor

Dynamic Scan, Easy Operation

LBY-XC40B/LBY-XC20B

Auto Dynamic Erythrocyte Sedimentation Rate (ESR) Monitor



Technical specifications

Testing Channel	40/20
Testing Time	60min/30min
Scanning Interval	2.5min
Accuracy	±5%
Channel Consistency	5%
Repeatability Error	≤3%
Report Unit	mm/h

Features

- Dynamic scanning of erythrocyte sedimentation process
- Automatically calibration to the test result at 25°C
- Auto induction of sample loading and auto testing
- Extremely easy operation: load and testing

PRECIL
A Mindray Company

Beijing Precil Instrument Co., Ltd.

Building No5, Shangdi Qunying Technology and Science Park, Haidian District, Beijing, P.R.China(100085)
TEL: 0086-10-62971818 62987758 FAX: 0086-10-62987761 <http://www.precil.com.cn>